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EN LEVENHUK 850B/870T/D870T BIOLOGICAL MICROSCOPES

General Information

Levenhuk 800 Series biological microscopes are safe for health, life and property of the consumer and the environment when properly used, and meets the requirements of international standards. These microscopes are designed for observing transparent objects in the transmitted light using the bright field method. They perform well in clinical research and tests; teaching demonstrations; bacterioscopy and cytoscopy in medical and health establishments, laboratories, universities, and schools; and may be used for scientific research in agriculture and microbiology.

The kit includes

- microscope body
- eyepiece head: binocular (Levenhuk 850B), trinocular (Levenhuk 870T and D870T)
- plan achromatic objectives: 4x, 10x, 40xs, 100xs (oil)
- eyepieces: WF10x/22mm (2pcs)
- blue filter
- vial of immersion oil
- dust cover
- 110–220V power cable
- 3W LED
- user manual and warranty

The Levenhuk D870T kit additionally includes: 8Mpx digital camera, a USB cable and a CD with software.

Optional accessories

You can increase the possibilities of your microscope by equipping it with the optional accessories: eyepieces of different magnification, eyepieces with reticle, dark-field condenser, LCD microscope camera, etc. Visit our website to choose and buy optional accessories: www.levenhuk.com

Microscope assembly

1. Eyepiece head. Loosen the locking screw on the bottom of the microscope body. Install the eyepiece head and tighten the locking screw.
2. Eyepiece. There is a special ring to fix the eyepiece firmly in the tube of eyepiece head. Make sure to install the proper end of the eyepiece in the tube.
3. Installing objective lenses. To install the objective lenses, screw them in by turning clockwise. Attach objective lenses to the revolving nosepiece in ascending order of magnification.
4. Condenser. The condenser was mounted, focused and centered at the factory. You can adjust the condenser, if needed, by using the centering screws.
5. Digital camera (Levenhuk D870T). Insert the camera into the C-mount and tighten the locking screw. Then, insert the C-mount into the trinocular tube.

Usage

- Illumination. Connect your microscope to a power outlet. Turn the power on; the power indicator will come on. You can change the level of brightness with the brightness adjustment wheel.
- Iris diaphragm is incorporated into the condenser. Adjust the iris diaphragm to match the numerical aperture of different objectives. Open or close the iris diaphragm with the adjusting knob. When the knob is on the right, the iris diaphragm is fully opened.
- Placing the specimen. Place the specimen on the stage and cover it with glass. Make sure to fix the specimen securely with the stage dials. Move the stage in horizontal and vertical directions to point the objective lens at the specimen.
- Focusing. Select the objective lens with the lowest magnification rotating the revolving nosepiece. Raise the stage turning the coarse focus knob. While observing through the objective lens, turn the coarse focus knob slowly to let the stage down until you see an image of the observed specimen. Turning the fine focus knob, make the image sharp. CAUTION: The objective should not touch the specimen; otherwise, the objective and/or the specimen might be damaged.
- Adjusting the eyepiece head. Adjust the eyepiece tubes to make the view in the left and right eyepiece even. Start the diopter adjustment with "0" value and continue until the image becomes sharp.
- Immersion oil. For observing at high magnifications with a 100x objective lens, you have to use immersion oil. Place a drop of immersion oil on the specimen. Turn the revolving nosepiece to place the 100x objective lens in its working position. Raise the microscope stage with a coarse focus knob until the bottom of the objective lens touches the oil drop. Move the revolving nosepiece back and forth to get rid of air bubbles in the oil. Then, fix the objective lens in its working position. Make sure that oil fills the entire space between the objective lens and the specimen.
- Changing the bulb. The microscope is equipped with built-in LED illumination. Change the bulb only if necessary.

Digital camera

Levenhuk D870T comes with a 8Mpx digital camera. The camera allows you to observe specimens in fine detail and true colors on your PC monitor and save images on the hard drive. Connect the camera to your PC via a USB-cable that also acts as a power cable.

Microscope Specifications

	850B	870T	D870T
Magnification, x	40–1000 (1600/2000)*		
Tube length	Infinity (∞)		
Observation head	Binocular head	Trinocular head	
	Eyepiece diameter: 30mm, Gemel binocular head (Siedentopf, rotatable 360°), eye relief: 40mm, 180° of rotating capability (at a 65mm interpupillary distance), 30° inclined, interpupillary distance: 48–75mm (180° of rotating capability when the interpupillary distance is 65mm, microscope head magnification: 1x		
Camera	–	–	+
Eyepieces, x/field, mm	10x/22mm; diopter adjustment: ±5 diopters on each eyepiece		
Nosepiece	4 objectives		
Optical design	Infinity achromatic objectives, parfocal distance: 45mm		
Objectives, x/aperture	4x/0.1; 10x/0.25; 40x/0.65; 100x/1.25		
Mechanical stage	Two-axis mechanical stage, no X-axis positioning rack, stage size: 216mm×150mm, moving range: 75mm×55mm		
Condenser	Abbe condenser, N.A. 1.25, center-adjustable, with adjustable aperture diaphragm, height-adjustable, dovetail mount		
Focusing mechanism	Coaxial coarse & fine focusing knobs (left side & both sides), coarse focusing tension adjusting knob, coarse focusing lock knob		
Light source	3W LED		
AC power supply			
voltage, V	220±22		
frequency, Hz	50		
Operating temperature range, °C	+5... +35		
Operating humidity range, %	20...80		

* The magnification of the microscope can be increased by using additional (optional) eyepieces and objectives. The manufacturer reserves the right to make changes to the product range and specifications without prior notice.

Camera Specifications

Megapixels	8
Video recording	+
Mounting location	eyepiece tube (replaces the eyepiece)
Image format	*.bmp, *.jpg, *.png, *.tif
Exposure	ERS
White balance	auto/manual
Exposure control	auto/manual
Software	Levenhuk

Software image control	image size, brightness, exposure control
Port	USB 2.0, 480Mb
System requirements	Windows XP/Vista/7/8/10/11 (32 and 64 bit), Mac OS 10.12, Linux Ubuntu 14.04
Camera power supply	USB cable

Care and maintenance

Never, under any circumstances, look directly at the Sun, another bright source of light or at a laser through this device, as this may cause PERMANENT RETINAL DAMAGE and may lead to BLINDNESS. Take necessary precautions when using the device with children or others who have not read or do not fully understand these instructions. After unpacking your microscope and before using it for the first time, check for the integrity and durability of every component and connection. Do not try to disassemble the device on your own for any reason. For repairs and cleaning of any kind, please contact your local specialized service center. Protect the device from sudden impact and excessive mechanical force. Do not apply excessive pressure when adjusting focus. Do not overtighten the locking screws. Do not touch the optical surfaces with your fingers. To clean the device exterior, use only special Levenhuk cleaning wipes and optics cleaning tools. Do not use any corrosive or acetone-based fluids to clean the optics. Abrasive particles, such as sand, should not be wiped off the lenses, but instead blown off or brushed away with a soft brush. Do not use the device for lengthy periods of time or leave it unattended in direct sunlight. Keep the device away from water and high humidity environments. Be careful during your observations; always place the dust cover to protect the device from dust and stains when not in use or when being stored. If do not use your microscope for a long time, store the objective lenses and eyepieces separately from the microscope. Store the device in a dry, cool place away from hazardous acids and other chemicals, away from heaters, open fires, and other sources of high temperatures. When using the microscope, try not to use it near flammable materials or substances (benzene, paper, cardboard, plastic, etc.), as the base may become hot during operation and cause a fire. Always disconnect the microscope from a power source before opening the base or changing the bulb. Regardless of the bulb type (halogen or incandescent), give it time to cool down before attempting to replace it, and always replace it with the same type of a bulb. Always use the power supply with the proper voltage, i.e., as specified in the specifications of a microscope. Connecting the instrument to a different power outlet may damage the electric circuitry of the microscope, burn out the bulb, or even cause a short circuit. Children should use the device under adult supervision only. Seek medical advice immediately if a small part or a battery is swallowed.

Levenhuk International Lifetime Warranty

All Levenhuk telescopes, microscopes, binoculars, and other optical products, except for their accessories, carry a lifetime warranty against defects in materials and workmanship. A lifetime warranty is a guarantee on the lifetime of the product on the market. All Levenhuk accessories are warranted to be free of defects in materials and workmanship for six months from the date of purchase. The warranty entitles you to free repair or replacement of your Levenhuk product in any country where a Levenhuk representative office is located, provided that all terms and conditions of the warranty are met.

For further details, please visit: www.levenhuk.com/warranty

If warranty problems arise or if you need assistance in using the product, please contact your local Levenhuk branch.

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Levenhuk® is a registered trademarks of Levenhuk, Inc.
© 2006–2023 Levenhuk, Inc. All rights reserved.
Levenhuk Inc. (USA): 928 E 124th Ave. Ste D, Tampa, FL 33612, USA, +1-813-468-3001, contact_us@levenhuk.com